



Toileting Plus

The Oxford® Toileting Plus is an easy to fit sling designed specifically to facilitate the toileting process, providing excellent access to clothing. Made from a heavy-duty polyester material, it offers both excellent support and durability. It is not a general-purpose sling and will only suit around 25% of clients.

- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Adjustable, elasticated body belt for added support and security
- ✓ Padded under-arm and leg sections for added comfort
- ✓ Discreet colour coded tab on top section to denote size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed
(From seated position only)



Chair to Commode/Toilet

Size	Part No.	SWL
XS	SL1751	227 kg
S	SL1752	227 kg
M	SL1753	227 kg
L	SL1754	227 kg
XL	SL1755	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.