

High Back

The Oxford® High Back sling is designed for lifting and transferring clients from either a seated or recumbent position. Made from heavy-duty polyester and with integrated head support, it offers good support and reassurance for the client.

- ✓ Removable 'bones' to the back section offer added support
- ✓ Padded leg sections for increased comfort
- ✓ Positioning pocket enables easier application
- ✓ Discreet colour coded tab on top section to denote size
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SL1763	227 kg
S	SL1764	227 kg
M	SL1765	227 kg
L	SL1766	227 kg
XL	SL1767	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.