


MultiFit SL

The Oxford® MultiFit SL is a development of the standard MultiFit sling. It incorporates over-sized, deluxe 'paddle' style leg supports for added comfort and support, particularly for clients with larger thighs. Integrated low-friction rip-stop nylon to the leg and spine sections ensures easy application and helps to reduce the risks associated with tissue viability.



- ✓ Optional detachable head support
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size		Part No.	SWL
1	XS	SL1540	227 kg
2		SL1541	227 kg
3	S	SL1542	227 kg
4		SL1543	227 kg
5	M	SL1544	227 kg
6		SL1545	227 kg
7		SL1546	227 kg
8		SL1547	227 kg
9	L	SL1548	227 kg
10		SL1549	227 kg
11	XL	SL1550	284 kg
12		SL1551	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.