


MultiFit Reflex

The Oxford® MultiFit Reflex sling incorporates 'reflex' shoulder straps that provide automatic adjustment for clients who have a lack of seat control or are prone to extensor spasm or mild to moderate involuntary movement. Integrated low-friction rip-stop nylon to the leg and spine sections ensures easy application and helps reduce the risks associated with tissue viability.

- ✓ Integrated head support for added support
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size		Part No.	SWL
1	XS	SL1701	227 kg
2		SL1702	227 kg
3	S	SL1703	227 kg
4		SL1704	227 kg
5	M	SL1705	227 kg
6		SL1706	227 kg
7		SL1707	227 kg
8		SL1708	227 kg
9	L	SL1709	227 kg
10		SL1710	227 kg
11	XL	SL1711	284 kg
12		SL1712	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.