





MultiFit SL Reflex

The Oxford® MultiFit SL Reflex sling incorporates 'reflex' shoulder straps that provide automatic adjustment for clients who have a lack of seat control or are prone to extensor spasm. It incorporates over-sized, deluxe 'paddle' style leg supports for added comfort and support, particularly for clients with larger thighs. Integrated low-friction rip-stop nylon to the leg and spine sections ensures easy application and helps reduce the risks associated with tissue viability.



- ✓ Integrated head support for added support
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



Size		Part No.	SWL
1	XS	SL1721	227 kg
2		SL1722	227 kg
3	S	SL1723	227 kg
4		SL1724	227 kg
5	M	SL1725	227 kg
6		SL1726	227 kg
7		SL1727	227 kg
8		SL1728	227 kg
9	L	SL1729	227 kg
10		SL1730	227 kg
11	XL	SL1731	284 kg
12		SL1732	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.