

Quickfit Spacer


The Oxford® Quickfit Spacer sling is an easy to fit, general purpose straight leg sling, designed to suit the majority of clients. Made from a cushioned and breathable 3D woven 'Spacer' material with high air permeability, offering high levels of comfort.



- ✓ Padded leg sections for increased comfort
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SPECIAL	227 kg
S	SL1401	227 kg
M	SL1402	227 kg
L	SL1403	227 kg
XL	SPECIAL	227 kg

