


# UltraFine Straight Leg

The Oxford® UltraFine Straight Leg sling is made from a cushioned and breathable 3D woven 'Spacer' material with high air permeability and pressure redistribution properties. Integrated, low-friction rip-stop nylon to the leg sections ensures easier application.

- ✓ Available with or without head support
- ✓ Detachable leg and shoulder straps
- ✓ Colour coded binding tab on top section to denote size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



Size	Part No. (Without HS)	Part No. (With HS)	SWL
XS	SL1561	SL1821	227 kg
S	SL1562	SL1822	227 kg
M	SL1563	SL1823	227 kg
L	SL1564	SL1824	227 kg
XL	SL1565	SL1825	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.