


UltraFine Deluxe Leg

The Oxford® UltraFine Deluxe Leg sling incorporates flared leg sections for increased patient comfort and support. Made from a cushioned and breathable 3D woven 'Spacer' material with high air permeability, offering high levels of comfort. Integrated, low-friction rip-stop nylon to the leg sections ensures easier application.

- ✓ Available with or without head support
- ✓ Suitable for amputees following careful assessment
- ✓ Detachable leg and shoulder straps
- ✓ Colour coded binding tab on top section to denote size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



WITH HEAD SUPPORT




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No. (Without HS)	Part No. (With HS)	SWL
XS	SL1571	SL1827	227 kg
S	SL1572	SL1828	227 kg
M	SL1573	SL1829	227 kg
L	SL1574	SL1830	227 kg
XL	SL1575	SL1831	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.