


UltraFine Reflex

The Oxford® UltraFine Reflex sling is made from a cushioned and breathable 3D woven ‘Spacer’ material with high air permeability, offering high levels of comfort. With counter-balance straps to the shoulder sections, this sling is particularly helpful where a patient is prone to moderate involuntary movement or extension. Integrated, low-friction rip-stop nylon to the leg sections ensures easier application.



- ✓ Head support and body belt provides added support and security
- ✓ Suitable for amputees following careful assessment
- ✓ Colour coded binding tab on top section to denote size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SL1530	227 kg
S	SL1531	227 kg
M	SL1532	227 kg
L	SL1533	227 kg
XL	SL1534	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.