

Quickfit Deluxe

The Quickfit Deluxe sling is a development of the standard Quickfit design and is suitable for the vast majority of patients. Contoured leg sections provide increased hip and thigh support and it can be used for amputees following careful assessment.

- ✓ Available in polyester solid or net material options
- ✓ Available with or without head support
- ✓ Optional padded leg sections for increased comfort
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Can be customised to fulfil patient specific needs



WITH HEAD SUPPORT



Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No. (Solid)	Part No. (Net)	SWL
XS	SPECIAL	SPECIAL	227 kg
S	SL1050	SL1051	227 kg
M	SL1053	SL1054	227 kg
L	SL1056	SL1057	227 kg
XL	SPECIAL	SPECIAL	227 kg

Size	Part No. (Padded)	Part No. (Padded HS)	SWL
XS	SPECIAL	SPECIAL	227 kg
S	SL1050P	SL1060P	227 kg
M	SL1053P	SL1063P	227 kg
L	SL1056P	SL1066P	227 kg
XL	SPECIAL	SPECIAL	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.