


UniFit

The Oxford® UniFit is an easy to fit, general-purpose sling with divided leg sections, suitable for the vast majority of patients. Padded leg sections as standard provide excellent comfort and support for the client.


- ✓ Available in polyester or net material options as standard
- ✓ Twin leg strap attachment points may offer increased comfort
- ✓ Positioning handles allow easy patient repositioning
- ✓ Colour coded binding indicates sling size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No. (Polyester)	Part No. (Net)	SWL
XS	SL1630	SL1841	227 kg
S	SL1631	SL1842	227 kg
M	SL1632	SL1843	227 kg
L	SL1633	SL1844	227 kg
XL	SL1634	SL1845	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.