

# UniFit Deluxe

The Oxford® UniFit Deluxe is an easy to fit, general-purpose sling suitable for most basic transfer requirements. The contoured leg sections, complete with padding and hip tapes as standard, offer excellent comfort and reduced pressure in the upper leg area.

- ✓ Available in polyester solid or net material options
- ✓ Suitable for amputees following careful assessment
- ✓ Positioning handles allow easy patient repositioning
- ✓ Colour coded binding indicates sling size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair  
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No. (Polyester)	Part No. (Net)	SWL
XS	SL1791	SL1796	227 kg
S	SL1792	SL1797	227 kg
M	SL1793	SL1798	227 kg
L	SL1794	SL1799	227 kg
XL	SL1795	SL1800	284 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.