


Comfort Access Padded

The Oxford® Comfort Access Padded sling, with head support, has been designed specifically to facilitate the toileting process, allowing easy access to clothing. For use in conjunction with a 4-point positioning cradle, the patented Securi3 clip system ensures no inadvertent detachment of the sling. It is not a general purpose sling and will only suit approximately 25% of clients.

- ✓ Fully padded with integrated head support for additional support
- ✓ Adjustable body belt for added support and security
- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow patient repositioning
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed
(From seated position only)



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SPECIAL	227 kg
S	SL25576	227 kg
M	SL25577	227 kg
L	SL25578	227 kg
XL	SPECIAL	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.