

Toileting Plus

The Oxford® Toileting Plus is an easy to fit sling designed specifically to facilitate the toileting process, providing excellent access to clothing. Made from a heavy-duty polyester material, it offers both excellent support and durability. It is not a general-purpose sling and will only suit around 25% of clients.

- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Adjustable, elasticated body belt for added support and security
- ✓ Padded under-arm and leg sections for added comfort
- ✓ Discreet colour coded tab on top section to denote size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed
(From seated position only)



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

| Size | Part No. | SWL |
|------|----------|--------|
| XS | SL1751 | 227 kg |
| S | SL1752 | 227 kg |
| M | SL1753 | 227 kg |
| L | SL1754 | 227 kg |
| XL | SL1755 | 284 kg |



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.