


# Comfort Disposable

The Oxford® Comfort Disposable sling is an easy fit, contoured sling designed to accommodate 85-90% of clients. For use in conjunction with a 4-point positioning cradle, the patented Securi3 clip system ensures no inadvertent detachment of the sling. In the interests of avoiding cross-infection, this sling is specifically design for single patient use only.

- ✓ Patented Securi3 clip attachment points
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Modesty loop helps protect patient dignity during transfer




Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SL25610	227 kg
S	SL25611	227 kg
M	SL25612	227 kg
L	SL25613	227 kg
XL	SL25614	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.