

UltraFine Lite Straight Leg

The Oxford® UltraFine Lite Straight Leg is a divided leg sling, made from a super-soft microfibre polyester, offering excellent levels of comfort. Integrated, low-friction rip-stop nylon to the leg sections ensures easy application.


- ✓ Available with or without head support
- ✓ Detachable leg and shoulder straps
- ✓ Colour coded binding tab on top section to denote size
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



WITH HEAD SUPPORT



Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair
(Head Support recommended)



Rehab/Assisted Standing

Size	Part No. (Without HS)	Part No. (With HS)	SWL
XS	SL1886	SL1891	227 kg
S	SL1887	SL1892	227 kg
M	SL1888	SL1893	227 kg
L	SL1889	SL1894	227 kg
XL	SL1890	SL1895	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.