

# Patient Handling Belt


The Patient Handling Belt offers carers a means of providing standing assistance to patients. It can be used independently or in combination with a compatible manual standing/transfer device, e.g. the Oxford Switch. It may also be used to help with walking rehabilitation/gait training activity following careful assessment.



- ✓ Side handles for standing assistance and support
- ✓ Adjustable belt enables use with the Oxford Switch
- ✓ Non-slip inner layer ensures secure fit
- ✓ Elasticated top band for increased comfort
- ✓ Colour coded binding indicates belt size
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed  
*(From/to seated position)*



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

| Size | Part No. | SWL    |
|------|----------|--------|
| XS   | SPECIAL  | 180 kg |
| S    | SL3151   | 180 kg |
| M    | SL3152   | 180 kg |
| L    | SL3153   | 180 kg |
| XL   | SPECIAL  | 180 kg |



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.