



Core Sling Dimensions: Oxford Access Padded

Access Padded			
Sling Size	Dimension A	Dimension B	Dimension C
XS	600mm	260mm	570mm
S	700mm	280mm	660mm
M	850mm	310mm	890mm
L	1050mm	330mm	1000mm
XL	1300mm	390mm	1160mm



*All dimensions are subject to manufacturing tolerances.

Access Padded with Head Support (Inc. Comfort/Clip Models)				
Sling Size	Dimension A	Dimension B	Dimension C	Dimension D
XS	590mm	270mm	1180mm	440mm
S	690mm	300mm	1360mm	490mm
M	830mm	310mm	1550mm	540mm
L	1060mm	280mm	1600mm	590mm
XL	1260mm	340mm	1860mm	670mm



*All dimensions are subject to manufacturing tolerances.

Custom Sling Service:

Joerns Healthcare offer a custom sling service should you have any non-standard or patient specific requirements. Please contact Joerns Healthcare for more information and guidance. Extended lead times apply.