

## Core Sling Dimensions: Oxford Patient Handling Belt

| Patient Handling Belt |             |             |
|-----------------------|-------------|-------------|
| Sling Size            | Dimension A | Dimension B |
| XS                    | -           | -           |
| S                     | 1210mm      | 250mm       |
| M                     | 1320mm      | 250mm       |
| L                     | 1480mm      | 250mm       |
| XL                    | -           | -           |



\*All dimensions are subject to manufacturing tolerances.

### Custom Sling Service:

Joerns Healthcare offer a custom sling service should you have any non-standard or patient specific requirements. Please contact Joerns Healthcare for more information and guidance. Extended lead times apply.