



Deluxe Standing with Clips

The Oxford® Deluxe Standing sling with Clips is an easy fit sling suitable for standing and supporting. It is not a general purpose sling and clients should be carefully assessed for suitability. It will not suit double amputees or lifting from a horizontal position and clients must have a good degree of weight bearing ability.


- ✓ Fully padded design offers excellent comfort
- ✓ Patented Securi3 clip attachment points
- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Adjustable waist strap for added support and security
- ✓ Can be customised to fulfil patient specific needs




Chair to Chair



Bed to Chair/Chair to Bed
(From/to seated position)



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair



Rehab/Assisted Standing

Size	Part No.	SWL
XS	SPECIAL	200 kg
S	SL3162	200 kg
M	SL3163	200 kg
L	SL3164	200 kg
XL	SPECIAL	200 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.